

Table 3. Formulation Compositions of Solid Microemulsion Preconcentrate (SM) and Solid Dispersion (SD)

Component	Solid Microemulsion Preconcentrate		PEG-Surfactant-based Solid Dispersion	
Compound 1 (% w/w)	4	8	4	8
Capmul PG8 (% w/w)	28.8	27.6	—	—
Cremophor EL (% w/w)	28.8	27.6	—	—
Polysorbate 80 (% w/w)	—	—	24	23
PEG 3350 (% w/w)	38.4	36.8	72	69

other parts. The SEM pictures in Figure 6 show that the liquid preconcentrate exists as a separate phase trapped into solid PEG 3350 clusters. Confocal Raman microscopy described in the next section showed a similar result (Fig. 6c and d): two phases were detected with crystalline needle-like and irregular shaped domains of PEG 3350 (visualized in green) and a liquid phase (visualized in blue) consisting of Capmul PG8 and Cremophor EL located between PEG 3350 domains. In separate experiments, hot-stage light microscopy studies (not described in this report) showed that the clusters and nodes disappeared upon heating to 60°C, and reappeared upon cooling below 50°C.

Characterization of Drug-Loaded Solid Microemulsion Preconcentrates by Confocal Raman Microscopy

The compositions of the solid microemulsion preconcentrate with 4% and 8% drug loads selected for further evaluation as potential dosage forms for **I** are given in Table 3. For comparison, compositions of prototype solid dispersions developed previously⁶ using PEG 3350-polysorbate 80 mixtures are also given in the table.

One well-established technique for identification and characterization of the solid state of pharmaceutical specimens is Raman spectroscopy. It is a robust and reliable method yielding information of high quality with a low limit of detection. In addition it requires minimal sample preparation and allows for high throughput screening and automation.²⁷

In combination with an optical microscope, Raman spectroscopy can be implemented to map a specimen and to provide spatially resolved chemical information on the underlying species.^{28–32} In recent years, progress in laser equipment, development of high quality holographic filters and gratings as well as advances in

charge couple device (CCD) technology has led to an advent of Raman microscopy. Employing the principle of confocality, modern Raman microspectrometers allow rapid chemical mapping with high spatial and spectral resolution. It is possible to reveal the exact composition and spatial distribution of complex mixtures of drug substances and excipients at a resolution in the submicrometer range, at very low concentrations and at a pixel resolution of up to several hundred thousand pixels.^{31,33}

In the present investigation, Raman spectra of the components of solid microemulsions in their pure form were first collected (Fig. 7) and used as the basis for the deconvolution of the experimental data cubes using classical least squares




Figure 7. Reference Raman spectra obtained from pure substances. These spectra were employed as basis for the deconvolution of the hyperspectral data cubes using classical least squares. From top to bottom: drug **I** (red), PEG3350 (green), Capmul PG8 (blue), Cremophor EL (black), Polysorbate 80 (orange).